

Energy Losses in Bends

Instruction Manual

F1-22

ISSUE 8

November 2012

038214

Table of Contents

Copyright and Trademarks	1
General Overview	2
Equipment Diagrams.....	3
Important Safety Information.....	4
Introduction.....	4
Electrical Safety.....	4
Description	5
Overview.....	5
Installation	6
Installing the Equipment and Commissioning.....	6
Operation	8
Operating the Equipment.....	8
Equipment Specifications.....	9
Environmental Conditions.....	9
Routine Maintenance	10
Responsibility	10
General.....	10
Laboratory Teaching Exercises.....	11
Index to Exercises	11
Nomenclature	11
Technical Data.....	12
Exercise A.....	13
Operating the Optional Software (F1-301).....	17
Contact Details for Further Information.....	24

Disclaimer

This document and all the information contained within it is proprietary to Armfield Limited. This document must not be used for any purpose other than that for which it is supplied and its contents must not be reproduced, modified, adapted, published, translated or disclosed to any third party, in whole or in part, without the prior written permission of Armfield Limited.

Should you have any queries or comments, please contact the Armfield Customer Support helpdesk (Monday to Thursday: 0830 – 1730 and Friday: 0830 - 1300 UK time). Contact details are as follows:

United Kingdom	International
(0) 1425 478781 (calls charged at local rate)	+44 (0) 1425 478781 (international rates apply)
Email: support@armfield.co.uk	
Fax: +44 (0) 1425 470916	

Copyright and Trademarks

Copyright © 2012 Armfield Limited. All rights reserved.

Any technical documentation made available by Armfield Limited is the copyright work of Armfield Limited and wholly owned by Armfield Limited.

Brands and product names mentioned in this manual may be trademarks or registered trademarks of their respective companies and are hereby acknowledged.

General Overview

Fluid mechanics has developed as an analytical discipline from the application of the classical laws of statics, dynamics and thermodynamics, to situations in which fluids can be treated as continuous media. The particular laws involved are those of the conservation of mass, energy and momentum and, in each application, these laws may be simplified in an attempt to describe quantitatively the behaviour of the fluid.

The hydraulics bench service module, F1-10, provides the necessary facilities to support a comprehensive range of hydraulic models each of which is designed to demonstrate a particular aspect of hydraulic theory.

The specific hydraulic model that we are concerned with for this experiment is the Energy Losses in Bends and Fittings Rig, F1-22. A full description of the apparatus is given later in these texts.

Equipment Diagrams




Figure 1: F1-22 Energy Losses in Bends Apparatus

Important Safety Information

Introduction

Before proceeding to operate the equipment described in this text we wish to alert you to potential hazards so that they may be avoided.

Although designed for safe operation, any laboratory equipment may involve processes or procedures which are potentially hazardous. The major potential hazards associated with this particular equipment are listed below.

- Injury through misuse
- Injury from electric shock
- Damage to clothing
- Risk of infection due to lack of cleanliness

Accidents can be avoided provided that equipment is **regularly maintained** and **staff** and **students** are made aware of potential hazards list of general safety rules is included in the F1 Product Manual to assist **staff** and **students** in this regard. The list is not intended to be fully comprehensive but for guidance only.

Please refer to the notes in the F1 Product Manual regarding the Control of Substances Hazardous to Health Regulations.

Electrical Safety

The F1-10 Service Bench operates from a mains voltage electrical supply. The equipment is designed and manufactured in accordance with appropriate regulations relating to the use of electricity. Similarly, it is assumed that regulations applying to the operation of electrical equipment are observed by the end user.

However, to give increased operator protection, Armfield Ltd have incorporated a Residual Current Device (RCD, alternatively called an Earth Leakage Circuit Breaker or ELCB) as an integral part of the service bench. If through misuse or accident the equipment becomes electrically dangerous, an RCD will switch off the electrical supply and reduce the severity of any electric shock received by an operator to a level which, under normal circumstances, will not cause injury to that person.

Check that the RCD is operating correctly by pressing the TEST button. The circuit breaker **MUST** trip when the button is pressed. Failure to trip means that the operator is not protected and the equipment must be checked and repaired by a competent electrician before it is used.

Description

Where necessary, refer to the drawings in the Equipment Diagrams section.

Overview

The accessory is designed to be positioned on the side channels of the hydraulics bench top channel.

The following fittings are connected in a series configuration to allow direct comparison:

- long bend
- area enlargement
- area contraction
- elbow bend
- short bend
- valve fitting
- mitre bend

Flow rate through the circuit is controlled by a flow control valve.

Pressure tapplings in the circuit are connected to a twelve bank manometer, which incorporates an air inlet/outlet valve in the top manifold. An air bleed screw facilitates connection to a hand pump. This enables the levels in the manometer bank to be adjusted to a convenient level to suit the system static pressure.

A clamp which closes off the tapplings to the mitre bend is introduced when experiments on the valve fitting are required. A differential pressure gauge gives a direct reading of losses through the gate valve.

Installation

Installing the Equipment and Commissioning

The Energy Losses in Bends apparatus is supplied ready for use and only requires connection to the F1-10 Hydraulics Bench as described below.

Carefully remove the components from the cardboard packaging. Retain the packaging for future use.

Locate the apparatus over the moulded channel in the top of the bench and ensure that the base plate is horizontal. Adjust the feet on the base plate if necessary.

Connect the flexible inlet tube at the left hand end to the quick release fitting in the bed of the channel.

Place the free end of the flexible outlet tube in the volumetric tank of the bench.

Fully open the gate valve and the outlet flow control valve at the right hand end of the apparatus.

Close the bench flow control valve then start the service pump.

Gradually open the bench flow control valve and allow the pipework to fill with water until all air has been expelled from the pipework.

In order to bleed air from the pressure tapping points and the manometers close both the bench flow control valve and the outlet flow control valve and open the air bleed screw. Remove the cap from the adjacent air inlet/outlet connection. Connect a length of small bore tubing from the air valve to the volumetric tank. Now, open the bench flow control valve and allow flow through the manometers to purge all air from them; then, tighten the air bleed screw and partly open both the bench valve and the outlet flow control valve. Next, open the air bleed screw slightly to allow air to enter the top of the manometers, re-tighten the screw when the manometer levels reach a mid height.

Gradually increase the volume flowrate until the pattern just fills the range of the manometer (adjust the bench flow control valve and the outlet flow control valve in combination to maintain all of the readings within the range of the manometer). If the pattern is too low on the manometer open the bench flow control valve to increase the static pressure. If the pattern is too high open the outlet flow control valve on the apparatus to reduce the static pressure.

These levels can be adjusted further by using the air bleed screw and the hand pump supplied. The air bleed screw controls the air flow through the air valve, so when using the hand pump, the bleed screw must be open. To retain the hand pump pressure in the system, the screw must be closed after pumping.

If the levels in the manometer are too high then the hand pump can be used to pressurise the top manifold. All levels will decrease simultaneously but retain the appropriate differentials.

If the levels are too low then the hand pump should be disconnected and the air bleed screw opened briefly to reduce the pressure in the top manifold. Alternatively

the outlet flow control valve can be closed to raise the static pressure in the system which will raise all levels simultaneously.

If the level in any manometer tube is allowed to drop too low then air will enter the bottom manifold. If the level in any manometer tube is too high then water will enter the top manifold and flow into adjacent tubes.

Note: If the static pressure in the system is excessive, Eg. with the bench flow control valve fully open and the outlet flow control valve almost closed, it will not be possible to use the hand pump to lower the levels in the manometer tubes. The valves should be adjusted to provide the required flowrate at a lower static pressure.

In operation the pressure drop across each fitting is compared with the volume flowrate which is measured using the volumetric measuring tank and a stopwatch (not supplied).

To check the operation of the differential pressure gauge associated with the gate valve, close off the flexible connecting tubes to the mitre bend pressure tapings using the clamps supplied before closing the gate valve (to prevent air being drawn into the system).

Open the bench flow control valve and the outlet flow control valve. As the gate valve is closed the differential pressure across the valve will be displayed on the gauge.

Close the bench flow control valve then switch off the service pump.

The F1-22 Energy Losses in Bends apparatus is ready for use.

Operation

Operating the Equipment

See Laboratory Teaching Exercises for details on operating the equipment.

Equipment Specifications

Environmental Conditions

This equipment has been designed for operation in the following environmental conditions. Operation outside of these conditions may result reduced performance, damage to the equipment or hazard to the operator.

- a. Indoor use;
- b. Altitude up to 2000m;
- c. Temperature 5°C to 40°C;
- d. Maximum relative humidity 80% for temperatures up to 31°C, decreasing linearly to 50% relative humidity at 40°C;
- e. Mains supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage;
- f. Transient over-voltages typically present on the MAINS supply;

Note: The normal level of transient over-voltages is impulse withstand (over-voltage) category II of IEC 60364-4-443;

- g. Pollution degree 2.

Normally only nonconductive pollution occurs.

Temporary conductivity caused by condensation is to be expected.

Typical of an office or laboratory environment.

Routine Maintenance

Responsibility

To preserve the life and efficient operation of the equipment it is important that the equipment is properly maintained. Regular maintenance of the equipment is the responsibility of the end user and must be performed by qualified personnel who understand the operation of the equipment.

General

Little maintenance is required but it is important to drain all water from the pipework when not in use. The air tapping on the top manifold of the manometer can be used to eject all water from the manometer tubes by opening the bleed valve then using the hand pump if necessary.

Any manometer tube which does not fill with water or is slow to fill or empty indicates that the tapping in the pipework or the connection at the base of the manometer tube is blocked or partially blocked. Disconnect the flexible connecting tube between the pipe fitting and the manometer. Blowing through the tapping will usually dislodge any foreign body.

A few drops of wetting agent introduced into the manometer tubes will reduce the meniscus with the glass wall and improve measurement accuracy.

Laboratory Teaching Exercises

Index to Exercises

Exercise A

Nomenclature

Name	Unit	Symbol	Type	Definition
Diameter of Test Pipe	m	d	Given	Internal diameter of test pipe.
Volume Collected	m ³	V	Measured	Taken from scale on hydraulics bench
Time to Collect	s	t	Measured	Time taken to collect the known volume of water in the hydraulics bench.
Temp of Water	°C		Measured	Temperature of water under test.
Kinematic Viscosity	(m ² /s)	ν	Measured	Taken from the table Kinematic Viscosity of Water at Atmospheric Pressure
Manometer	m	h_1	Measured	Measured value from the appropriate manometer. Note: Scale calibrated in mm.
Manometer	m	h_2	Measured	Measured value from the appropriate manometer. Note: Scale calibrated in mm.
Timed Flow Rate	m ³ /s	Q_t	Calculated	$Q_t = \frac{V}{t} = \frac{\text{Volume Collected}}{\text{Time to Collect}}$
Velocity	m/s	v	Calculated	Velocity of fluid in duct = $\frac{Q_t}{A} = \frac{\text{Flow Rate}}{\text{Area}}$
Dynamic Head	m		Calculated	$\frac{v^2}{2g}$
Reynolds Number		Re	Calculated	$Re = \frac{vd}{\nu}$

Head Loss/ Pressure Drop	m	Δh	Calculated	Head loss across the fitting $\Delta h = h_2 - h_1$
Loss Coefficient		K	Calculated	$K = \frac{\Delta h}{\frac{v^2}{2g}}$ See Theory

Technical Data

The following dimensions from the equipment are used in the appropriate calculations. If required these values may be checked as part of the experimental procedure and replaced with your own measurements.

Internal diameter of pipework: $d = 0.0183\text{m}$

Internal diameter of pipework at enlargement outlet and contraction inlet:

$$d = 0.0240\text{m}$$

Exercise A

Objective

To determine the loss factors for flow through a range of pipe fittings including bends, a contraction, an enlargement and a gate-valve.

Method

By measurement of head differences across each of a number of fittings connected in series, over a range of steady flows.

Equipment

In order to complete the demonstration we need a number of pieces of equipment.

- The F1-10 Hydraulics Bench which allows us to measure flow by timed volume collection.
- The F1-22 Energy Losses in Bends and Fittings Apparatus.
- A stopwatch to allow us to determine the flow rate of water (not supplied).
- Clamps for pressure tapping connection tubes.
- Spirit level (not supplied).
- Thermometer (not supplied).

Theory

The energy loss which occurs in a pipe fitting (so-called secondary loss) is commonly expressed in terms of a head loss (h , metres) in the form:

$$\Delta h = \frac{Kv^2}{2g}$$

where K = the loss coefficient and v = mean velocity of flow into the fitting.

Because of the complexity of flow in many fittings, K is usually determined by experiment. For the pipe fitting experiment, the head loss is calculated from two manometer readings, taken before and after each fitting, and K is then determined as

$$K = \frac{\Delta h}{v^2 / 2g}$$

Due to the change in pipe cross-sectional area through the enlargement and contraction, the system experiences an additional change in static pressure. This change can be calculated as

$$\frac{v_1^2}{2g} - \frac{v_2^2}{2g}$$

To eliminate the effects of this area change on the measured head losses, this value should be added to the head loss readings for the enlargement and the contraction.

Note that $(h_1 - h_2)$ will be negative for the enlargement and negative for the contraction.

$$\frac{v_1^2}{2g} - \frac{v_2^2}{2g} \text{ will be}$$

For the gate valve experiment, pressure difference before and after the gate is measured directly using a pressure gauge. This can then be converted to an equivalent head loss using the equation

$$1 \text{ bar} = 10.2 \text{m water}$$

The loss coefficient may then be calculated as above for the gate valve.

Reynolds number is a dimensionless number used to compare flow characteristics. A full investigation of Reynolds number, and typical flow variation as it changes, can be found in experiment F1-20 (Reynolds' Apparatus) or in a suitable fluids text book.

Equipment Set Up

Set up the losses apparatus on the hydraulic bench so that its base is horizontal (this is necessary for accurate height measurements from the manometers). Connect the test rig inlet to the bench flow supply and run the outlet extension tube to the volumetric tank and secure it in place.

Open the bench valve, the gate valve and the flow control valve and start the pump to fill the test rig with water. In order to bleed air from pressure tapping points and the manometers close both the bench valve and the test rig flow control valve and open the air bleed screw and remove the cap from the adjacent air valve. Connect a length of small bore tubing from the air valve to the volumetric tank. Now, open the bench valve and allow flow through the manometers to purge all air from them; then, tighten the air bleed screw and partly open both the bench valve and the test rig flow control valve. Next, open the air bleed screw slightly to allow air to enter the top of the manometers, re-tighten the screw when the manometer levels reach a convenient height.

Check that all manometer levels are on scale at the maximum volume flow rate required (approximately 17 litres/minute). These levels can be adjusted further by using the air bleed screw and the hand pump supplied. The air bleed screw controls the air flow through the air valve, so when using the hand pump, the bleed screw must be open. To retain the hand pump pressure in the system, the screw must be closed after pumping.

Results

It is not possible to make measurements on all fittings simultaneously and, therefore, it is necessary to run two separate tests.

Test 1 measures losses across all pipe fittings except the gate valve, which should be kept fully open. Adjust the flow from the bench control valve and, at a given flow rate, take height readings from all of the manometers after the levels have steadied. In order to determine the volume flow rate, you should carry out a timed volume collection using the volumetric tank. This is achieved by closing the ball valve and measuring (with a stopwatch) time taken to accumulate a known volume of fluid in the tank, which is read from the sight glass. You should collect fluid for at least one minute to minimise timing errors.

Repeat this procedure to give a total of at least five sets of measurements over a flow range from approximately 8 - 17 litres per minute. Measure the outflow water temperature at the lowest flow rate; this together with the table detailing the Kinematic Viscosity of Water at Atmospheric Pressure is used to determine the Reynolds number.

Test 2 measures losses across the gate valve only. Clamp off the connecting tubes to the mitre bend pressure tapplings (to prevent air being drawn into the system). Start with the gate valve closed and open fully both the bench valve and the test rig flow control valve. Now open the gate valve by approximately 50% of one turn (after taking up any backlash). For each of at least 5 flow rates, measure pressure drop across the valve from the pressure gauge; adjust the flow rate by use of the test rig flow control valve. Once measurements have started, do not adjust the gate valve. Determine the volume flow rate by timed collection.

Repeat this procedure for the gate valve opened by approximately 70% of one turn and then approximately 80% of one turn.

FITTING	Manometer h_1 m	Manometer h_2 m	Head Loss $h_1 - h_2$ m	Vol V m^3/s	Time t s	Flow Rate Q_t m^3/s	Velocity v m/s	$\frac{2}{25}$	K
MITRE									
ELBOW									
SHORT BEND									
ENLARGEMENT									
CONTRACTION									
GATE VALVE	Gauge Reading =								

Conclusion

For Test 1, plot graphs of head loss (Δh) against dynamic head, and K against volume flow rate Q_t .

For Test 2, plot graphs of equivalent head loss (Δh) against dynamic head, and K against Q_t .

Comment on any relationships noticed. What is the dependence of head losses across pipe fittings upon velocity?

Examining the Reynolds number obtained, are the flows laminar or turbulent?

Is it justifiable to treat the loss coefficient as constant for a given fitting?

In Test 2, how does the loss coefficient for a gate valve vary with the extent of opening the valve?

Kinematic Viscosity of Water at Atmospheric Pressure


Temperature (degrees C)	Kinematic Viscosity ν ($10^{-6} \times \text{m}^2/\text{s}$)	Temperature (degrees C)	Kinematic Viscosity ν ($10^{-6} \times \text{m}^2/\text{s}$)
0	1.793	25	0.893
1	1.732	26	0.873
2	1.674	27	0.854
3	1.619	28	0.836
4	1.568	29	0.818
5	1.520	30	0.802
6	1.474	31	0.785
7	1.429	32	0.769
8	1.386	33	0.753
9	1.346	34	0.738
10	1.307	35	0.724
11	1.270	36	0.711
12	1.235	37	0.697
13	1.201	38	0.684
14	1.169	39	0.671
15	1.138	40	0.658
16	1.108	45	0.602
17	1.080	50	0.554
18	1.053	55	0.511
19	1.027	60	0.476
20	1.002	65	0.443
21	0.978	70	0.413
22	0.955	75	0.386
23	0.933	80	0.363
24	0.911	85	0.342

Eg. At 20°C the kinematic viscosity of water is $1.002 \times 10^{-6} \text{m}^2/\text{s}$.


Operating the Optional Software (F1-301)

Note: The diagrams in this section are included as typical examples and may not relate specifically to an individual product.

The Armfield Software is a powerful Educational and Data Logging tool with a wide range of features. Some of the major features are highlighted below, to assist users, but full details on the software and how to use it are provided in the presentations and Help text incorporated in the Software. Help on Using the Software or Using the Equipment is available by clicking the appropriate topic in the **Help** drop-down menu from the upper toolbar when operating the software as shown:



Load the software. If multiple experiments are available then a menu will be displayed listing the options. Wait for the presentation screen to open fully as shown:



Presentation Screen - Basics and Navigation


As stated above, the software starts with the Presentation Screen displayed. The user is met by a simple presentation which gives them an overview of the capabilities of the equipment and software and explains in simple terms how to navigate around the software and summarizes the major facilities complete with direct links to detailed context sensitive 'help' texts.

To view the presentations click **Next** or click the required topic in the left hand pane as appropriate. Click **More** while displaying any of the topics to display a Help index related to that topic.

To return to the Presentation screen at any time click the View Presentation icon



from the main tool bar or click **Presentation** from the dropdown menu as shown:



For more detailed information about the presentations refer to the **Help** available via the upper toolbar when operating the software.

Toolbar

A toolbar is displayed at the top of the screen at all times, so users can jump immediately to the facility they require, as shown:




The upper menu expands as a dropdown menu when the cursor is placed over a name.

The lower row of icons (standard for all Armfield Software) allows a particular function to be selected. To aid recognition, pop-up text names appear when the cursor is placed over the icon.


Mimic Diagram

The Mimic Diagram is the most commonly used screen and gives a pictorial representation of the equipment, with boxes to enter measurements from the equipment, display any calculated variables etc. directly in engineering units.

To view the Mimic Diagram click the View Diagram icon  from the main tool bar or click **Diagram** from the **View** drop-down menu as shown:



A Mimic diagram is displayed, similar to the diagram as shown:



The details in the diagram will vary depending on the equipment chosen if multiple experiments are available.


Manual data input boxes with a coloured background allow measured variables, constants such as Orifice Cd and Atmospheric Pressure, as appropriate, to be changed by over-typing the default value. After typing the value press the Return key or click on a different box to enter the value.

In addition to measured variables such as Volume, Time, Temperature or Pressure, calculated data such as Discharge / Volume flowrate, Headloss etc are continuously displayed in data boxes with a white background. These are automatically updated and cannot be changed by the user.

After entering a complete set of data from measurements on the equipment click on

the  icon to save the set of results before entering another set.

The mimic diagram associated with some products includes the facility to select different experiments or different accessories, usually on the left hand side of the screen, as shown:



Clicking on the appropriate accessory or exercise will change the associated mimic diagram, table, graphs etc to suit the exercise being performed.


Data Logging Facilities in the Mimic Diagram

Armfield software designed for manual entry of measured variables does not include automatic data logging facilities and these options are greyed out where not

appropriate. When manually entering data the  icon simply saves the set of entered data into a spreadsheet as described above.

Tabular Display

To view the Table screen click the View Table icon  from the main tool bar or click **Table** from the View dropdown menu as shown:



The data is displayed in a tabular format, similar to the screen as shown:

Sample Number	Rate	Accumulated Pressure	Temperature	Conversion	Rate	Rate	Rate	Rate	Rate	Rate	Rate
1	...	...	...	...	...	...	...	...	...	...	...
2	...	...	...	...	...	...	...	...	...	...	...
3	...	...	...	...	...	...	...	...	...	...	...
4	...	...	...	...	...	...	...	...	...	...	...
5	...	...	...	...	...	...	...	...	...	...	...
6	...	...	...	...	...	...	...	...	...	...	...
7	...	...	...	...	...	...	...	...	...	...	...
8	...	...	...	...	...	...	...	...	...	...	...
9	...	...	...	...	...	...	...	...	...	...	...
10	...	...	...	...	...	...	...	...	...	...	...
11	...	...	...	...	...	...	...	...	...	...	...
12	...	...	...	...	...	...	...	...	...	...	...

As the data is sampled, it is stored in spreadsheet format, updated each time the data is sampled. The table also contains columns for the calculated values.

New sheets can be added to the spreadsheet for different data runs by clicking the




icon from the main toolbar. Sheets can be renamed by double clicking on the sheet name at the bottom left corner of the screen (initially Run 1, Run 2 etc) then entering the required name.

For more detailed information about Data Logging and changing the settings within the software refer to the **Help** available via the upper toolbar when operating the software.


Graphical Display

When several samples have been recorded, they can be viewed in graphical format.

To view the data in Graphical format click the View graph icon  from the main tool bar or click **Graph** from the **View** drop-down menu as shown:




The results are displayed in a graphical format as shown:




(The actual graph displayed will depend on the product selected and the exercise that is being conducted, the data that has been logged and the parameter(s) that has been selected).

Powerful and flexible graph plotting tools are available in the software, allowing the user full choice over what is displayed, including dual y axes, points or lines, displaying data from different runs, etc. Formatting and scaling is done automatically by default, but can be changed manually if required.

To change the data displayed on the Graph click **Graph Data** from the **Format** dropdown menu as shown:



The available parameters (Series of data) are displayed in the left hand pane as shown:




Two axes are available for plotting, allowing series with different scaling to be presented on the same x axis.

To select a series for plotting, click the appropriate series in the left pane so that it is highlighted then click the appropriate right-facing arrow to move the series into one of the windows in the right hand pane. Multiple series with the same scaling can be plotted simultaneously by moving them all into the same window in the right pane.

To remove a series from the graph, click the appropriate series in the right pane so that it is highlighted then click the appropriate left-facing arrow to move the series into the left pane.

The X-Axis Content is chosen by default to suit the exercise. The content can be changed if appropriate by opening the drop down menu at the top of the window.

The format of the graphs, scaling of the axes etc. can be changed if required by clicking **Graph** in the **Format** drop-down menu as shown:



For more detailed information about changing these settings refer to the **Help** available via the upper toolbar when operating the software.

Contact Details for Further Information

Main Office: Armfield Limited

Bridge House
West Street
Ringwood
Hampshire
England BH24 1DY

Tel: +44 (0)1425 478781

Fax: +44 (0)1425 470916

Email: sales@armfield.co.uk
support@armfield.co.uk

Web: <http://www.armfield.co.uk>

US Office: Armfield Inc.

9 Trenton - Lakewood Road
Clarksburg, NJ 08510

Tel/Fax: (609) 208 2800

Email: info@armfieldinc.com